


Install Anaconda on facility's station


1. Go to "sfgclients" on Desktop
2. Go to: SCF > DHSI folder, and click on Anaconda installer package



3. Follow the instruction.
4. One thing you need to be careful about is the installation destination:
Choose "Install for me only" for Destination



5. Now Anaconda is installed on your account: home directory (the house icon in Finder window). You can launch Anaconda-Navigator from your_Netlink home directory/anaconda/Anaconda-Navigator.



Note: At the time of logoff, everything on your home directory will be erased. This means that “anaconda” and “AnacondaProject” will be gone when you logoff. You may want to save the project on your **USB** or “Desktop” or “My Document” volumes (icons) which you see on your local Desktop.